



219 North Arkansas Street
Rogers, AR 72756
Phone: (479) 636-4361; 800-643-3426
Fax: (479) 636-3562
Email: biosales@pel-freez.com
Website: www.pelfreez-bio.com

PRODUCT DATA SHEET

RABBIT BLOOD, NON-STERILE, WITH SODIUM HEPARIN

Product Codes:

31132 -1

Pack size:

500 mL

Description:

Rabbit whole blood is collected from individual, healthy, young rabbits (New Zealand White) 8 - 12 weeks of age, mixed sex at time of slaughter, into Sodium Heparin (7500 units/ml) at a ratio of 500:1, bottled and stored on wet ice prior to shipment.

Physical State:

Fresh, non-frozen liquid

Testing:

No analytical testing is performed on this product. The blood is visually inspected at the time of bottling to ensure the absence of clots.

Packaging, shipping/storage:

Packaging

Polypropylene bottle with lid

Storage Temperature

2 – 8 °C

Shipping Conditions

This product is shipped fresh on wet ice at 2 – 8 °C.

Expiration

This product is perishable. The exact shelf life will depend on the user's application.