



219 North Arkansas Street
Rogers, AR 72756
Phone: (479) 636-4361; 800-643-3426
Fax: (479) 636-3562
Email: biosales@pel-freez.com
Website: www.pelfreez-bio.com

PRODUCT DATA SHEET

PIG GAMMA GLOBULIN

Product Codes:

29315-1

Pack size:

1 g

Description:

Porcine gamma globulin is purified from normal porcine serum or plasma using a combination of existing methods developed in our laboratory.

Physical State:

Off-white, fine, flaky powder

Testing:

Serum Protein Electrophoresis: $\geq 95\%$

Protein by Weight: O.D. 280 nm $E^{1\%1.39}$: $\geq 85\%$ - Residual salts, lipids, etc., comprise the rest of the weight.

Moisture Content: FIO

Packaging, shipping/storage:

Packaging

Polypropylene bottle

Storage Temperature

2 – 8 °C

Shipping Conditions

Ambient temperature

Expiration

Product quality is guaranteed to meet Pel-Freez Biologicals' specifications for 1 year from the date of receipt by the customer as long as the product is stored in accordance with the indicated storage conditions.

Application Notes: Reportedly useful as a stabilizer/blocker for EIA and RIA systems, IHC assays, immunoprecipitation, Western Blot and ELISA experiments, and protein standards.