

219 North Arkansas Street Rogers, AR 72756

Phone: (479) 636-4361; 800-643-3426 Fax: (479) 636-3562

> Email: biosales@pel-freez.com Website: www.pelfreez-bio.com

## PRODUCT DATA SHEET

## RABBIT SERUM, STERILE, NON-HEMOLYZED

 Product Codes:
 Pack size:

 31125-1
 100 mL

 31125-5
 500 mL

### **Description:**

Whole blood is collected from young healthy rabbits (New Zealand White or Californian, 8 – 12 weeks of age, either sex and fasted) at the time of slaughter, refrigerated within 6 hours, processed within 24 hours and frozen at -20 °C. The serum is thawed, pooled, passed through a 0.22 micron filter, bottled, and returned to -20 °C.

### **Physical State:**

Frozen liquid

#### Testing:

Total Protein: 5.0 - 7.0 g/dL Hemolysis: ≤ 20 mg/dL

pH: 7.0 - 8.5

Sterility: Passed, no growth at 14 days

## Packaging, shipping/storage:

# **Packaging**

PETG bottle

### **Storage Temperature**

-10 °C or below

### **Shipping Conditions**

Product is shipped in frozen state, either with dry ice or in a refrigerated cargo carrier.

## **Expiration**

Product quality is guaranteed to meet Pel-Freez Biologicals' specifications for 1 year from the date of receipt by the customer as long as the product is stored in accordance with the indicated storage conditions.

<u>Application Notes</u>: Reportedly useful in IHC assays, protein purification, and as cell culture and growth media.

#### References:

Chromosomal Map Locatin of the Apha-Hemolysin Structural Gene in *Staphylococcus aureus* NCTC 8325: P.A. Pattee, Dept of Microbiology and Dept of Genetics, Iowa State University, Ames, Iowa 50011: Infection and Immunicty, Nov. 1986 Vol. 54, No.2, p. 593-59.

P/N: 74030 Rev 01 Page 1 of 1