



219 North Arkansas Street  
Rogers, AR 72756  
Phone: (479) 636-4361; 800-643-3426  
Fax: (479) 636-3562  
Email: [biosales@pel-freez.com](mailto:biosales@pel-freez.com)  
Website: [www.pelfreez-bio.com](http://www.pelfreez-bio.com)

## PRODUCT DATA SHEET

### PIG SERUM, NON-STERILE, TRACE-HEMOLYZED

**Product Codes:**

39519-5

**Pack size:**

500 mL

**Description:**

Whole blood is collected from healthy pigs (mixed breed, either sex) at the time of slaughter, refrigerated within 6 hours, processed within 24 hours and frozen at -20 °C.

**Physical State:**

Frozen liquid

**Testing:**

Hemolysis:  $\geq 20 \leq 50$  mg/dL

**Packaging, shipping/storage:**

**Packaging:**

Polypropylene bottle

**Storage Temperature:**

-20 °C or below

**Shipping Conditions:**

Dry ice

**Expiration:**

Product quality is guaranteed to meet Pel-Freez Biologicals' specifications for 1 year from the date of receipt by the customer as long as the product is stored in accordance with the indicated storage conditions.

**Application Notes:**

Reportedly useful in IHC assays and protein purification.