



219 North Arkansas Street
Rogers, AR 72756
Phone: (479) 636-4361; 800-643-3426
Fax: (479) 636-3562
Email: biosales@pel-freez.com
Website: www.pelfreez-bio.com

PRODUCT DATA SHEET

BOVINE SERUM, STERILE, NON-HEMOLYZED

Product Codes:

37225-5

Pack size:

500 mL

Description:

Bovine whole blood is collected from healthy animals (any breed, either sex, ≥ 1 yr old) refrigerated within 6 hours of collection, processed within 24 hours and frozen at -20°C . The serum is thawed and passed through a 0.22 micron filter, bottled and returned to -20°C .

Physical State:

Frozen liquid

Testing:

Hemolysis: ≤ 20 mg/dL

Total Protein: 6.3-8.9 g/dL

pH: 7.0-8.5

Sterility: Passed, no growth at 14 days

Packaging, shipping/storage:

Packaging:

Polyethylene glycol bottle

Storage Temperature:

-20°C or below

Shipping Conditions:

Dry ice

Expiration:

Product quality is guaranteed to meet Pel-Freez Biologicals' specifications for 1 year from the date of receipt by the customer as long as the product is stored in accordance with the indicated storage conditions.

Application Notes:

Reportedly useful in IHC assays, protein purification, and as a nutrient source in cell culture/growth media