



219 North Arkansas Street
Rogers, AR 72756
Phone: (479) 636-4361; 800-643-3426
Fax: (479) 636-3562
Email: biosales@pel-freez.com
Website: www.pelfreez-bio.com

PRODUCT DATA SHEET

BOVINE PLASMA, NON-STERILE, SODIUM HEPARIN

Product Codes

37142-1

Pack size:

500 mL

Description:

Bovine whole blood is collected using Sodium Heparin, 10 units/mL. The plasma is separated from the blood cellular components by centrifugation and bottled. The animals are ≥ 1 year old and come from the continuous 48 United States. The product was obtained from USDA inspected facilities.

Physical State:

Frozen liquid

Packaging, shipping/storage:

Packaging

Polypropylene bottles

Storage Temperature

-20 °C or below

Shipping Conditions

Dry ice

Expiration

Product quality is guaranteed to meet Pel-Freez Biologicals' specifications for 1 year from the date of receipt by the customer as long as the product is stored in accordance with the indicated storage conditions.

Application Notes:

Reportedly useful for protein purification