



219 North Arkansas Street
Rogers, AR 72756
Phone: (479) 636-4361; 800-643-3426
Fax: (479) 636-3562
Email: biosales@pel-freez.com
Website: www.pelfreez-bio.com

PRODUCT DATA SHEET

CHICKEN WHOLE BLOOD, NON-STERILE, SODIUM CITRATE

Product Codes

33130-1

Pack size:

500 mL

Description:

Healthy animals 6-8 weeks of age, mixed gender and breed; collected on wet ice following USDA guidelines in 3.8% Sodium Citrate solution at 1:10 dilution.

Physical State:

Fresh, non-frozen liquid

Testing:

No analytical testing is performed on this product. The blood is visually inspected at the time of bottling to ensure the absence of clots.

Packaging, shipping/storage:

Packaging

Polypropylene bottle

Storage Temperature

2 – 8 °C

Shipping Conditions

Shipped on day of collection on wet ice

Expiration:

This product is perishable. The exact shelf life will depend upon the user's application.