



219 North Arkansas Street
Rogers, AR 72756
Phone: (479) 636-4361; 800-643-3426
Fax: (479) 636-3562
Email: biosales@pel-freez.com
Website: www.pelfreez-bio.com

PRODUCT DATA SHEET

BOVINE WHOLE BLOOD, NON-STERILE, SODIUM EDTA

Product Codes

37131-1

Pack size:

500 mL

Description:

Healthy animals aged > 1yr, mixed gender and breed; collected on wet ice following USDA guidelines, using 2.5% Sodium Ethylenediaminetetraacetic Acid solution at 1:10 dilution

Physical State:

Fresh, non-frozen liquid

Testing:

No analytical testing is performed on this product. The blood is visually inspected at the time of bottling to ensure the absence of clots.

Packaging, shipping/storage:

Packaging

Polypropylene bottle

Storage Temperature

2 – 8 °C

Shipping Conditions

Shipped on the day of collection on wet ice.

Expiration:

This product is perishable. The exact shelf life will depend upon the user's application.