



219 North Arkansas Street  
Rogers, AR 72756  
Phone: (479) 636-4361; 800-643-3426  
Fax: (479) 636-3562  
Email: [biosales@pel-freez.com](mailto:biosales@pel-freez.com)  
Website: [www.pelfreez-bio.com](http://www.pelfreez-bio.com)

## PRODUCT DATA SHEET

### PIG WHOLE BLOOD, NON-STERILE, SODIUM EDTA

**Product Codes**

39531-1

**Pack size:**

500 mL

**Description:**

Pig blood, mixed age, gender, and breed is collected following USDA guidelines using 2.5% Sodium Ethylenediaminetetraacetic Acid solution at 1:10 dilution.

**Physical State:**

Fresh, non-frozen liquid

**Testing:**

No analytical testing is performed on this product. The blood is visually inspected at the time of bottling to ensure the absence of clots.

**Packaging, shipping/storage:**

**Packaging**

Polypropylene bottle

**Storage Temperature**

2 – 8 °C

**Shipping Conditions**

Shipped on the day of collection on wet ice.

**Expiration:**

This product is perishable. The exact shelf life will depend upon the user's application.